

Table 1 Primers with sequences used in this study

Primer Name	Sequence (5'-3')	Product Size
Cry2A-F	AGATTACCCCAGTTCCAGAT	500 bp
Cry2A-R	GTTCCCGAAGGACTTTCTAT	
GTG-F	CCCTGGTGACAAGTCCATCT	350 bp
GTG-R	CTGCACACCCATCTCTCTGA	
Cry1Ac-F	ACAGAAGACCCTTCAATATC	1 Kbp
Cry1Ac-R	GTTACCGAGTGAAGATGTAA	