

Table 1 Quantitative PCR procedures

Name of the gene	Pre denaturation	Denaturation	Annealing	Extension	Cycles
<i>GoMADS-box2</i>	94°C, 2 min	94°C, 30 s	58°C, 30 s	72°C, 30 s	35
<i>Goactin1</i>	94°C, 2 min	94°C, 30 s	60°C, 30 s	72°C, 30 s	35
<i>GoGAPDH1</i>	94°C, 2 min	94°C, 30 s	56°C, 30 s	72°C, 30 s	35